









GO SUPRA

| Nozzle name  | Powders  | Indications   |
|--|--|---|
| <br>AIR-N-GO® EASY SUPRA 120°<br><i>Supplied as standard</i>                | <br><b>"CLASSIC"</b><br>Sodium bicarbonate-based<br>Average particle size 76µm<br><i>FLAVOURS: Lemon, Peppermint, Raspberry, Cola and Neutral</i> | Finishing after scaling,<br>Preparation for whitening,<br>Cleaning of inter-proximal areas, fissures and crevices,<br>Preparation of the tooth surface before etching,<br>Cleaning of orthodontic brackets,<br>Removal of supra-gingival biofilm, plaque, stains and temporary cement residues. |
|  | <br><b>"PEARL"</b><br>Calcium carbonate-based<br>Average particle size 55µm<br><i>FLAVOUR: Neutral</i>  | "Ultra-fast" cleaning of difficult stains and composites with no risk of damaging them<br><b>Note: The absence of sodium makes it possible to treat sensitive patients or those following a salt-free diet</b>  |
| <br>AIR-N-GO® PERIO EASY<br><i>Supplied with the PERIO kit as an option</i> | <br><b>"PERIO"</b><br>Glycine-based<br>Average particle size 25µm<br><i>FLAVOUR: Neutral</i>  | Removal of the sub-gingival biofilm<br>Treatment of periodontal disease and peri-implantitis<br>Pockets 3mm - 8mm   |
|  |  | Treatment of periodontal disease and peri-implantitis<br>Pockets 8mm - 10mm   |
|  |  | Periodontal maintenance<br>Pockets up to 4mm  |
| <br>AIR-N-GO® PERIO MAINTENANCE<br><i>Supplied as an option</i>           |  |   |

For more information about good practice and the use of AIR-N-GO® easy, see our channel **You Tube**:



COMMERCIAL OFFER

Standard CI

**AIR-N-GO® easy case:** AIR-N-GO® easy handpiece with adaptor depending on the model (KaVo\*, Sirona\*, Bien-Air\*, W&H\*, Midwest\*), 1 x SUPRA 120° nozzle, 1x intro kit of 10 multi-flavour sticks of "CLASSIC" powder and 2 x sticks of "PEARL" powder, 1x AIR-N-GO® easy maintenance kit (syringe of silicone grease, applicator brushes for silicone grease, replacement o'rings, syringe and cannula, cleaning probes)

Optional

**PERIO easy Kit:** AIR-N-GO® PERIO easy nozzle and 1x bottle of 160 g AIR-N-GO® "PERIO" glycine powder  
**Supplied alone:** AIR-N-GO® PERIO nozzle  
**Supplied alone:** AIR-N-GO® PERIO Maintenance nozzle

\*KaVo, Bien-Air, W&H, Sirona and Midwest are registered trademarks.







Air polishing is becoming even simpler with

**AIR N GO**  
*easy*

Freshness  
Active control  
of periodontal  
disease

**ACTEON® EQUIPMENT** has created the new generation of portable air polishers, AIR-N-GO® easy with dual functions:

- "SUPRA" function" for supra-gingival polishing: daily prophylactic polishing.
- "PERIO" function" for sub-gingival polishing: maintenance, periodontal and peri-implant treatments.

This medical device allies simplicity with performance.

Bacterial decontamination  
Freedom of  
movement

Simplicity

Supra  
Perio



**A SINGLE HANDPIECE TO COVER  
A RANGE OF CARE AND TREATMENTS**

AIR-N-GO® easy combines essential functions with reliability to facilitate your everyday life:

- Fast start-up: direct connection to the high speed handpiece connector.
- Greater freedom of movement: 360° rotation, fast and accurate treatment of areas requiring care.
- Easier cleaning of the handpiece: main parts can be dismantled.

#### COMPLETE SOLUTION

4 nozzles for active control of the progress of periodontal disease.

Depending on your clinical needs, choose:

- Sodium bicarbonate-based "CLASSIC" powders
- Calcium carbonate-based "PEARL" powder
- Glycine-based "PERIO" powder.



Reliability  
Easy to convert  
Prevention

Flexibility

Comfort

Compatibility

Treatment of  
peri-implantitis

Power and  
performance

Patient safety



# airN GO nozzles:

## Your allies in preventing and fighting disease

These nozzles, which result from the most advanced fluid dynamics modelling techniques, provide continuous, accurate and controlled spray whatever the clinical context.

In your daily practice these AIR-N-GO® nozzles provide you with power and efficiency for the removal of supra- and sub-gingival biofilm.

### AIR-N-GO® SUPRA NOZZLE *(standard)*

- Tool suitable for any supra-gingival prophylactic treatment in combination with AIR-N-GO® "CLASSIC" or "PEARL" powders
- Guarantees a fast, efficient procedure for greater working comfort



*Clinical case Dr P. ROUAS (France)*



### AIR-N-GO® PERIO MAINTENANCE NOZZLE *(optional)*

- Reduced consumption of AIR-N-GO® "PERIO" fine powder particles in the maintenance phase
- Treatment of mucositis and return to a healthy clinical situation
- Action against the progression of periodontal disease



*Clinical case Dr C. LIM (U.K.)*



### AIR-N-GO® PERIO EASY NOZZLE *(optional PERIO easy kit)*

- Innovative design in a "shoehorn" shape providing easy sub-gingival insertion
- A real asset for non-surgical, sub-gingival treatment of teeth and implants (3 to 8 mm periodontal pockets)
- 5 seconds per face is all that is needed



*Clinical case Dr C. LIM (U.K.)*



### AIR-N-GO® PERIO NOZZLE *(optional)*

- Sub-gingival use for surgical and non-surgical treatments (tunnelling and open flaps)
- Powerful double lateral spray
- Decontamination of deep periodontal pockets (8-10 mm)



*Clinical case Dr A. PATEL (U.K.)*





# airN GO powders: customised solutions for your clinical needs

AIR-N-GO® powders have been specially designed and studied to produce an effective treatment that is gentle on all treated surfaces. These powders are compatible with the AIR-N-GO® easy handpiece and no change of instrument is required.

## Sodium bicarbonate-based AIR-N-GO® "CLASSIC" powders (≈ 76 µm)


- Complete range of supra-gingival powders which meet your patients' preferences
- Scrupulous cleaning giving instant results
- Action of removing coloured stains and bacterial plaque, even in areas that are difficult to access
- Better treatment acceptance: less bleeding and pain felt by the patient

### 100% NATURAL FRESHNESS




5 fresh flavours based on 100% natural aromas or essential oils:  
Neutral, Cola, Raspberry, Peppermint and Lemon


airN GO




Stratified particles with aggressive geometry (x1000)



Treated enamel shows numerous scratches (x300)



Particles with softened geometry (x1000)




Treated enamel: surface preserved without impacts or scratches (x300)

Source: ICMCB CNRS (Institute of Condensed Matter Chemistry - French National Centre for Scientific Research) S.E.M. photos

## Natural calcium carbonate-based AIR-N-GO® "PEARL" powder (≈ 55 µm)

- Ultra-active supra-gingival polishing and protection of fragile tissues
- Reduction in pH acidity caused by tissue inflammation; return to clinically healthy buccal flora
- Restoration of the natural shine and pearlyness of dental enamel
- Fast cleaning of difficult stains and composites in total safety

Even the most sensitive patients with advanced gingivitis will benefit from treatment provided with "PEARL".



AIR-N-GO® "PEARL" (x1000)

Source: ICMCB CNRS (Institute of Condensed Matter Chemistry - French National Centre for Scientific Research) S.E.M. photos

### HIGH PERFORMANCE




## Natural glycine-based AIR-N-GO® «PERIO» powder (≈ 25 µm)

- Exclusive use for sub-gingival treatments as a complement to the initial therapy
- Low abrasiveness preserving the treated surfaces (teeth and implants) and the most fragile anatomical areas
- Prevents the progression of the sub-gingival biofilm

### IN-DEPTH SUB-GINGIVAL THERAPY

Its unique properties with multiple advantages:

- Plaque removal: during periodontal maintenance on exposed root surfaces without sacrificing substance
- Reduction of the bacterial load in the periodontal pockets



AIR-N-GO® "PERIO" (x1000)

Source: ICMCB CNRS (Institute of Condensed Matter Chemistry - French National Centre for Scientific Research) S.E.M. photos